



**PRODUCT**  
CATALOGUE | 2018



# equipment biomedical services

USL IS MUCH MORE  
THAN JUST AN OFFICE AND A  
WAREHOUSE FULL OF STOCK.

IT IS A SERVICE THAT  
ENHANCES THE HEALTH OF  
THE COMMUNITY AROUND US.



## SEAMLESS SUPPORT FOR OUR CUSTOMERS.

Since 1984, USL has grown to be a significant supplier in the healthcare market of New Zealand. We supply a wide range of products to professionals nationwide in primary healthcare. We are suppliers to all 20 district health boards and public hospitals in new zealand. Today we have a significant presence in both consumer retail and industrial markets. USL Aesthetics is the newest addition to the USL Family, with an extensive line of products for the Aesthetically minded. By keeping with the USL philosophy we are known for, we are excited about this new direction and look forward to cultivating new partnerships and success in this market.

## OUR SERVICING DIVISION.

Our highly skilled Service Technicians are equipped to install, service and repair the full USL product range and other medical equipment. This includes lasers as well as other servicing and preventative maintenance requirements. USL Technicians work according to AS/NZS: 3551 2012 standard and our workshop complies with ISO9001. We are also the only authorised service agent for warranty repairs for all Welch Allyn products.

## IT'S ALL ABOUT PEOPLE.

USL has been built by creating seamless relationships. Since 1984 we've built enduring partnerships with customers, suppliers and our community. We set out to create the very best environment for our team so that they might deliver their unique skills to our markets.



## TABLE OF CONTENTS.

SCULPSURE .....	4 - 5
PICOSURE.....	6-7
AESTHETICS LASERS .....	8
PELLEVE .....	9
MOTIVA BREAST IMPLANTS .....	10-11
KELLER FUNNEL.....	12
LIPOELASTIC COMPRESSION GARMENTS.....	13
AESTHETICS CONSUMABLES.....	14-19
SIGVARIS & CONTACT DETAILS.....	20

# SculpSure™

25-minute  
**BODY  
CONTOURING**

## SCULPSURE


SculpSure is the first FDA-cleared laser treatment for non-invasive lipolysis of the flanks and abdomen. A state-of-the-art, controlled light-based technology, perfect for patients looking to get rid of stubborn fat, fast, effectively and in only 25mins. Over 90% patient satisfaction rate with 24% reduction in stubborn fat, and can be used in all skin types with no downtime.

### FLANK



### ABDOMEN





### TARGET

Selective wavelength to target fat cells below the dermis.

### DISRUPT

Highly efficient wavelength for controlled energy delivery.

### ELIMINATE

Over time, the body naturally eliminates the fat cells.

### CONTROLLED HYPERTHERMIC LASER TREATMENT

- Induces adipocyte injury by raising adipose temperature to a range of 42°C to 47°C
- No damage to surrounding tissue
- Comfortable and well tolerated procedure
- Proprietary energy modulation for efficacy and comfort

## LIGHT UP YOUR BUSINESS

Capitalize on the growing demand for energy-based devices like SculpSure.

**7 million** skin-tightening & body-shaping procedures were performed in 2014

**\$38 billion** earned from skin-tightening & body-shaping procedures performed in 2014

**15%** growth expected per year for the skin-tightening & body-shaping market through 2019



## FOCUS LENS ARRAY FOR SKIN REVITALISATION TREATMENTS

The FOCUS Lens Array harnesses the powerful photomechanical effect of the world's first and only picosecond aesthetic workstation. PicoSure FOCUS microscopically concentrates the PicoSure pulse to a precise depth, offering a skin revitalisation treatment on a wide range of skin types unlike any other.

The Focus Lens Array is ideal for skin revitalisation and pigmented lesions. This is accomplished by exposing less than 10% of the skin to areas of high fluence embedded in a low fluence background. PicoSure Focus treatments are ideal for patients seeking dramatic results with minimal downtime.





### The Hard Facts

- ▶ Powerful, dual-wavelength 755nm and 532nm picosecond laser system for the removal of tattoos of all colours
- ▶ Most comprehensive FDA cleared picosecond laser for tattoo removal, pigmentation, acne scarring, and wrinkles
- ▶ Largest compilation of clinically proven and peer reviewed papers by both local and international specialists
- ▶ Focus Lens Array - the only picosecond fractionated delivery system for skin revitalization (PicoToning™) with minimal to no patient downtime
- ▶ The only 755nm wavelength picosecond technology for the effective removal of blue and green inks
- ▶ Proven stability; proven technology; proven results! Asked for by name!



# Enjoy the Tattoos You Love...


Clear the Ones You Regret.  PicoSure








## PICOSURE

PicoSure, the world's first, fastest and effective picosecond aesthetic laser. This unparalleled breakthrough in laser technology delivers ultra-short pulse bursts of energy to the skin in trillionths of a second. Picosecond pulse width is 100 times shorter than nanosecond technology, enabling unmatched photomechanical impact for better clearance with fewer treatment, less fluence and without injury to the surrounding skin.





### Ink Response to Laser

Untreated Tissue/Ink	Q-Switched laser	PicoSure laser
 "Rocks"	 "Pebbles"	 "Sand"
		
	<b>Before</b>	<b>1 Month post 3 Tx</b>

Electron microscopy photos courtesy of H. R. Jallian, MD, Split tattoo image courtesy of Cynosure

## AESTHETIC LASERS

### MEDLITE C6 Q-SWITCHED ND:YAG

Gold-standard Medlite C6 features homogenous beam profile for uniform and consistent delivery of energy, and multiple wavelengths for efficient, reproducible treatments. Proven Elector Optic Q-Switched 532nm/1064nm technology enables consistent depth of penetration for fast effective results on all skin types for treatments including non-ablative skin resurfacing, melisma, Nevus of Ota, sun damage, tattoo removal, vascular lesion removal and more.



### EMERGE FRACTIONAL LASER

Introducing the lightweight, compact 1410nm Emerge Fractional Laser System that provides aesthetic treatments with little to no downtime and delivers deep, non-ablative laser micro-beams for fractional treatment. Ideal for treating fine lines and wrinkles, photo damage, general skin tone and texture and acne scarring.



### REVLITE SI

Globally considered the Gold-Standard in 1064nm/532nm Q-Switched Nd:YAG laser technology. Revlite features proprietary Photo Acoustic Technology Pulse™ to provide maximum power and versatility for a broad range of indications on all skin types such as non-ablative skin resurfacing, full colour spectrum tattoo removal, melasma, vascular lesion removal and more. The Smart Infinite™ hand piece and Photo Acoustic Technology Pulse provides enhanced power and versatility. True flat-top beam profile delivers high energy evenly over the skin's surface significantly reducing potential side effects.



### ICON AESTHETIC SYSTEM

Introducing the Icon Laser Light System, the next generation in laser light aesthetic systems with innovative and sophisticated functionality. The Icon allows you to offer unique FDA approved revenue generating treatment opportunities such as facial and leg veins, hair removal, skin resurfacing, wrinkle treatment, laser scar removal, stretch marks, skin rejuvenation and pigmentation.







## PELLEVE

The Pelleve System offers the proven benefits of radiofrequency (RF) to deliver excellent results to your patients. It is an excellent method of cutting and coagulating soft tissue, and FDA-cleared to treat facial wrinkles without surgery. With 4 waveforms (4 monopolar and 1 bipolar) and a wide selection of electrode tips, the user can easily configure the system for a wide array of applications. Non-invasive procedures include tightening skin above or below the eyes, wrinkles around the mouth, nasolabial folds crow's feet and more. Surgical procedures include biopsy, blepharoplasty, face lifts, lesion removal, telangiectasia, rhinophyma and keloids.



TREATMENT FOR MOLE REMOVAL



BEFORE



BEFORE

AFTER



4 MONTHS POST TREATMENT





No need to change your  
lifestyle with Motiva  
SilkSurface<sup>®</sup>, the surface of  
choice for active women

ZERO Double Capsules  
ZERO Late Seromas



## TRUETISSUE TECHNOLOGY - ERGONOMIX

Motiva Ergonomix™ are 100% filled implants that adapt shape after implantation to give the natural look of traditional “anatomical” implants, without the complications associated with rotation and implant hardness. TrueTissue Technology™ is the combination of a specific elastic elastomer shell and the special rheological properties of ProgressiveGel Ultima™ that combine with gravity to shift the maximum point of projection to the lower pole when the patient is in standing position. When the patient lies flat on her back, the implants will react in a similar way to a natural breast and the maximum point of projection will move closer to the middle point of her breast. Motiva Ergonomix™ implants adjust with gravity to the patient’s position, in order to provide a very natural result.

### ERGONOMIX


- We studied the rheological properties of a natural breast, in order to create the highly elastic Progressive Gel Ultima™ for a very soft, natural feeling, giving your patients the ultimate confidence. ProgressiveGel Ultima™ is made with the finest Nusil Technology silicone materials and does not fracture during implantation.

- TrueTissue Dynamics™ is the new gold standard for the natural behaviour of implants post-implantation. Both our TrueMonobloc® configuration and ProgressiveGel Ultima™ act in sync to mimic the natural behaviour of the breast. Giving patients the greatest of confidence when undertaking their daily lifestyle habits and activities.
- Our proprietary SilkSurface®, with its NanoSurface Technology™, does not adhere to surrounding tissues. Therefore, it allows the implant to adapt to the natural movement and eliminates the abrasion related to salt or sugar texturing, that can occur in physically active women.

### RANGE

- MOTIVA Implants are available in a wide selection of sizes for Round SilkSurface™ PLUS with Qid™ (Q-inside technology – the world’s first FDA cleared RFID micro-transponder for use in humans and the only one with a CE-Mark when used in a breast implant. Providing an electronic serial number from within the body for patient control and verification of their implant, for safety, well-being and in the event of a recall or adverse event), Round SilkSurface™ PLUS without Qid™ and Ergonomix Qid™

### ERGONOMIX IMPLANT





## EXPERIENCE THE EVOLUTION OF IMPLANT DELIVERY

### DISCOVER THE DIFFERENCE

Designed specifically to meet the needs of surgeons worldwide who employ textured, shaped and anatomical implants. The next generation KFTM2 delivers a clear alternative to traditional finger insertion techniques.



### BENEFITS

The Keller Funnel™2 is a cone-shaped, single-use device that provides a low friction insertion method for breast, calf, gluteal and pectoral silicone gel implantation procedures. Surgeons with experience using the funnel consistently report high satisfaction rates ([www.funnelstudy.com](http://www.funnelstudy.com)) and relate the funnel which allows them to:

- Shorten incision lengths
- Increase placement options
- Reduce force applied to the implant shell which may lessen rupture occurrence
- Employ a No-Touch technique believed by respondents to reduce infection rates
- Decrease operating room time

### 4 EASY STEPS

*The funnel is easy to use and requires knowing just 4 simple steps*



#### TRIM

Use guide to size funnel to implant



#### HYDRATE

Place funnel in sterile solution for a few seconds to hydrate the inner coating




#### LOAD

Place implant into funnel and confirm or adjust orientation



#### PROPEL

Squeeze behind implant and gently propel into surgical pocket




## DID YOU KNOW THAT A POST-OPERATIVE BRA IS ESSENTIAL FOR THE PERFECT RESULT OF AUGMENTATION?

### 1<sup>ST</sup> STEP

#### PI SPECIAL



- RECOMMENDED IMMEDIATELY AFTER AUGMENTATION
- RECOMMENDED USE: 2 WEEKS AFTER THE SURGERY



### 2<sup>ND</sup> STEP

#### PI UNIQUE



- RECOMMENDED AS A FOLLOW-UP CARE AFTER AUGMENTATION
- RECOMMENDED USE: 6 WEEKS FOLLOWING THE FIRST STEP



Invaluable aid for use during the first stage of post-surgery care to promote ideal breast shaping. The main advantage of this bra are the cups which are seamless, double-layered and heat-molded from elastic cotton material, which is very comfortable even during long-term wear. This design eliminates unwanted pressure or skin irritation while at the same time enabling easy manipulation. Specially reinforced cups ensure an even rate of compression and optimal breast support.

PI UNIQUE is the latest leading product in the area of compression bras from LIPOELASTIC. It is designed to be worn mainly by very demanding women who do not wish to compromise themselves even in the period of their convalescence. The shaped cups are designed for the second stage of post-surgery care. However, their design reveals much about their further potential use.

## SOFTFIL CANNULA

The SoftFil® Classic range of micro-cannulas is ideal for beginners as well as confirmed users. It has been developed after the renowned Soft Filling Technique, a standardised injection method accessible to all physicians.

- Broad series of Cannulas
- Adapted to every clinic application
- Designed for all fillers on the market
- Praised for it's quality and reliability
- Full clinical safety
- Minimal traumatism on tissues thanks to cannula's blunt tip



CODE	DESCRIPTION	UOM
SF2550L	Softfil Cannula 25g x 50mm	20 box
SF2740L	Softfil Cannula 27g x 40mm	20 box

## SALINE


Sodium Chloride injection.



CODE	DESCRIPTION	UOM
2417	Sodium Chloride for Injection 5ml	50 box
2388A	Sodium Chloride 0.9% Injection 20ml	50 box
2388B	Sodium Chloride 0.9% Injection 20ml	Each
2387	Sodium Chloride 0.9% Injection 10ml	50 box

## NEEDLES

Terumo hydermic needle.



## NUMBING CREAM

Emla numbing cream.



CODE	DESCRIPTION	UOM
4949	Emla Cream 5x5g	5 box

## CHLORHEXIDINE

Irrigation solution




CODE	DESCRIPTION	UOM
2420	Chlorhexidine 0.1% 1:2000 30ml	30 box
7460	PH Chorhex Acetate 0.05% 500ml Blue AHF7983	Each
7554	PH Chorhex Acetate 0.02% 100ml AHF7976	Each
2371	Chlorhexidine Aqueous 0.02% Irrigation 1:5000 500ml	Each
2372	Chlorhexidine 0.015% W/Cetrimide 0.15% Irrig 30ml	30 box
2373	Chlorhexidine 0.015% W/Cetrimide 0.15% Irrig 100ml	Each

CODE	DESCRIPTION	UOM
33322	Terumo Hypodermic Needle 18g x 1.5"	100 box
33225	Terumo Hypodermic Needle 19g x 1.5"	100 box
33227	Terumo Hypodermic Needle 20g x 1"	100 box
33229	Terumo Hypodermic Needle 20g x 1.5"	100 box
33230	Terumo Hypodermic Needle 21g x 1"	100 box
33232	Terumo Hypodermic Needle 21g x 1.5"	100 box
33234	Terumo Hypodermic Needle 22g x 1"	100 box
33236	Terumo Hypodermic Needle 22g x 1.5"	100 box
33235	Terumo Hypodermic Needle 22g x 2"	100 box
33237	Terumo Hypodermic Needle 23g x 1"	100 box
33238	Terumo Hypodermic Needle 23g x 1.25"	100 box
33239	Terumo Hypodermic Needle 24g x 1"	100 box
33240	Terumo Hypodermic Needle 25g x 5/8"	100 box
33241	Terumo Hypodermic Needle 25g x 1"	100 box
11887	Terumo Hypodermic Needle 26g x 0.5"	100 box
33246	Terumo Hypodermic Needle 27g x 0.5"	100 box

## SYRINGES

Terumo syringe.



CODE	DESCRIPTION	UOM
33012	Terumo Syringe 10ml L/L	100 box
33043	Terumo Syringe 10ml L/S	100 box
33040	Terumo Syringe 1ml Tuberculin	100 box
33013	Terumo Syringe 20ml L/L	100 box
33045	Terumo Syringe 20ml L/L Eccentric	100 box
33060	Terumo Syringe 30ml L/L	100 box
33010	Terumo Syringe 3ml L/L	100 box
33041	Terumo Syringe 3ml L/S	100 box
33011	Terumo Syringe 5ml L/L	100 box
33042	Terumo Syringe 5ml L/S	100 box

## ULTRASOUND GEL

AQUASONIC 100 is the most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required.



CODE	DESCRIPTION	UOM
3952	Ultrasound Gel 250ml	Each

## ASSORTED ANTIBACTERIAL OINTMENT

Betadine® Antiseptic Ointment is recommended for the treatment of common skin infections such as infected nails and tinea, as well as preventing infection in minor burns, cuts and abrasions. It is a bacterial, sporicidal, fungicidal and virucidal antiseptic.



CODE	DESCRIPTION	UOM
4812	Betadine Ointment 10% 25g	Each
10400	Bepanthen 50g	Each

## HEALTHE ANTIBACTERIAL OINTMENT

Recommended as a topical emollient for dry skin conditions such as eczema and dermatitis. This specially formulated dermatological cream makes it suitable as supportive therapy in conditions that have previously been treated with more potent corticosteroids. Non-perfumed and suitable for use on sensitive skin.



CODE	DESCRIPTION	UOM
10259	HE Fatty Cream	Each
10275	HE Antibacterial Hand Sanitiser Gel 375ml	Each
20001	HE Antibacterial Hand Sanitiser Gel 60ml	Each
10276	HE Chlorhexidine Hand Rub 500ml	Each
10763	HE Antiseptic Liquid Spray with Tea Tree Oil	Each
10351	HE Antiseptic Cream with Tea Tree Oil	Each

## USL NON-WOVEN SWABS

4 Ply Non Woven swabs in paper bag - Non Sterile.



CODE	DESCRIPTION	UOM
14390	USL Swab Non Woven Non Sterile 5x5cm 4ply	100 pkt
14391	USL Swab Non Woven Non Sterile 7.5x7.5cm 4ply	100 pkt
14392	USL Swab Non Woven Non Sterile 10x10cm 4ply	100 pkt

## USL DRESSING PACK

The USL wound dressing pack is an all in one pack ready to use for the cleaning and dressing of wounds. Able to be used as a base pack or a stand-alone pack.

- Easy to use
- Saves time, no need to collect the components individually
- Quality components at a low cost
- Individually wrapped and sterile

Contents:  
 1 x Tray 11x13.5x2.5cm (3 compartments)  
 6 x Non woven swabs 7.5 x 7.5cm  
 2 x Green forceps  
 1 x Yellow forcep  
 1 x Paper towel 42 x 40cm  
 1 x Sterile field 50 x 50cm



CODE	DESCRIPTION	UOM
14452	USL Dressing Pack #1	Each

## USL ASSORTED

Assorted USL branded product.



CODE	DESCRIPTION	UOM
14060	Apron Plastic LPDE White Disposable 71 x 117cm	100 box
14213	USL Applicator Cotton Single Tip 15cm	100 pk
14216	USL Tongue Depressor	100 box

## USL GAUZE SWABS

4 Ply Non Woven swabs in sterile pouch packaging 25 pouches per dispenser box - Sterile.



CODE	DESCRIPTION	UOM
14393	USL Swab Non Woven Sterile 5's 7.5x7.5cm Box 25	Box 25
14394	USL Swab Non Woven Sterile 2's 7.5x7.5cm Box 50	Box 50

## ALCOHOL SWABS

A skin cleansing swab saturated with 70% isopropyl alcohol.

- 2 ply swabs
- Box of 200



CODE	DESCRIPTION	UOM
12955	Alcohol Swabs	200 box

## USL CREPE BANDAGE

The cotton crepe bandage for support and retention. USL crepe bandage can be used as universal dressing retention bandages. All three bandaging techniques offer the benefits of support and relief of injuries while not inhibiting mobility.



CODE	DESCRIPTION	UOM
7355	USL Crepe Bandage 7.5cm x 4m	Each
7356	USL Crepe Bandage 5cm x 4m	Each
7357	USL Crepe Bandage 10cm x 4m	Each
7358	USL Crepe Bandage 15cm x 4m	Each



## USL NON-STERILE LATEX GLOVES

USL Examination Gloves are designed to meet and exceed international product standards using sophisticated technologies and innovations to provide uncompromising protection to its users.



CODE	DESCRIPTION	UOM
18884	USL Exam Glove Latex Powder Free - XS	100 box
18885	USL Exam Glove Latex Powder Free - S	100 box
18886	USL Exam Glove Latex Powder Free - M	100 box
18887	USL Exam Glove Latex Powder Free - L	100 box
18888	USL Exam Glove Latex Powder Free - XL	100 box

## USL UNI-GRIP TUBULAR BANDAGES

Comfortable, Effective Support Bandage for Sprains, Strains and Weak Joints.



CODE	DESCRIPTION	UOM
8985	USL Uni-Grip Tubular Bandage 10m Size A	Roll
8986	USL Uni-Grip Tubular Bandage 10m Size B	Roll
8987	USL Uni-Grip Tubular Bandage 10m Size C	Roll
8988	USL Uni-Grip Tubular Bandage 10m Size D	Roll
8989	USL Uni-Grip Tubular Bandage 10m Size E	Roll
8990	USL Uni-Grip Tubular Bandage 10m Size F	Roll
8991	USL Uni-Grip Tubular Bandage 10m Size G	Roll

## USL TRANSPARENT/PAPER TAPE

USL Paper Tape is a hypoallergenic adhesive paper tape. The tape is thin, flexible and easy to tear. Suitable for adhering first aid dressings.



## USL NON-STERILE NON-LATEX GLOVES

- Polymer Coated
- Powder Free
- Finger Textured Surface
- Single Use
- Bead Cuff
- Ambidextrous
- Chemotherapy Tested
- Colour: Blue
- Shelf Life: 3 Years



CODE	DESCRIPTION	UOM
19994K	USL Glove Nitrile Polymer Powder Free - XS	200 box
19995K	USL Glove Nitrile Polymer Powder Free - S	200 box
19996K	USL Glove Nitrile Polymer Powder Free - M	200 box
19997K	USL Glove Nitrile Polymer Powder Free - L	200 box
19998K	USL Glove Nitrile Polymer Powder Free - XL	200 box

## USL MEDICAL INSTANT COLD PACKS

- Instant and Rapid Cooling
- Reduces swelling and relieves heat stroke discomfort
- Size: 15cm x 22cm



CODE	DESCRIPTION	UOM
11070	USL Medical Instant Cold Pack 12cm x 15cm	Each
11071	USL Medical Instant Cold Pack 12cm x 15cm	12/Box
12885	USL Medical Instant Cold Pack 15cm x 22cm	Each
12885A	USL Medical Instant Cold Pack 15cm x 22cm	6/Box
13050	USL Medical Wrap 'N' Gel Re-Useable Cold Pack	Each
13050A	USL Medical Wrap 'N' Gel Re-Useable Cold Pack	6/Box

CODE	DESCRIPTION	UOM
9803	USL Paper Tape 12mm x 9.14m White	24 box
9799	USL Tape Transparent 12mm x 9.14m	24 box
9801	USL Tape Transparent 25mm x 9.14m	12 box
2041	Tape Durapore 12mm x 9.14m	Roll
2042	Tape Durapore 12mm x 9.14m	Roll
1290	Omniplast Rigid Medical Tape 2.5cm x 5m	Roll
1295	Omniplast Rigid Medical Tape 2.5cm x 9m	Roll
2021	Tape Transpore 12mm x 9.14m	Roll
2022	Tape Transpore 12mm x 9.14m	Roll

## LASTOTEL CONFORMING BANDAGE

The comfortable dressing retention bandage. Lastotel® is a dressing retention bandage with an extensibility of approximately 140%. It does not restrict or impede blood flow or joint movement when applied correctly.



CODE	DESCRIPTION	UOM
1322	Lastotel Cohesive Bandages 10cm x 4m Flesh	Each
7183	Lastotel Cohesive Bandages 12cm x 4m Flesh	Each

## PEHA-HAFT COHESIVE BANDAGE

A cohesive fixation bandage providing light support

- Peha-haft is composed of 67% viscose, 32% polyamide and 1% cohesive coating. It has a longitudinal extensibility of 80% and a transverse extensibility of 50%.



CODE	DESCRIPTION	UOM
1346	Peha-haft Elastic Bandage 2.5cm x 4m	8 box
1347	Peha-haft Elastic Bandage 4cm x 4m	Each
1348	Peha-haft Elastic Bandage 6cm x 4m	Each
1349	Peha-haft Elastic Bandage 8cm x 4m	Each
14523	Peha-haft Cohesive 4cm x 20m	Roll
1351	Peha-haft Elastic Bandage 6cm x 20m	Each
1352	Peha-haft Elastic Bandage 8cm x 20m	Each

## COBAN 2 COMPRESSION BANDAGE

Coban™ 2 Compression System materials were engineered with Intelligent Compression Dynamics to provide a conforming, inelastic sleeve with the required stiffness to support muscle pump to support lymph flow and reduce Lymphoedema.



CODE	DESCRIPTION	UOM
11370	Coban 2 Compression 5cm 20022	Roll
11371	Coban 2 Compression 10cm 20024	Roll
11732	Coban 2 Compression 15cm 20026	Roll
11362	Coban 2 Lite Compression 2.5cm 20721	Roll
11363	Coban 2 Lite Compression 7.5cm 20723	Roll
11364	Coban 2 Lite Compression 10cm 20724	Roll
11365	Coban 2 Lite Compression 15cm 20726	Roll
10721	Coban 2 Lite Complete System	Each

## PATIENT GOWN

- Unique wrap around design and modest length for optimal coverage
- Soft light weight synthetic fabric offering greater comfort
- Properly placed, easy to reach ties for patients to use without assistance
- Easy to process, quick drying
- Durable and stain resistant
- No fold system compliant



CODE	DESCRIPTION	UOM
90125	Patient Gown - Snowflake	Each

## PATIENT ISOLATION GOWN

USL Non-woven Isolation Gown with knitted cuffed sleeve.



CODE	DESCRIPTION	UOM
14678	Patient Isolation Gown - Yellow	Each

### LASTOPRESS ELASTIC COHESIVE BANDAGE

A lightweight permanent elastic cohesive bandage. This product is self adhering and is used for support, compression and securing applications.



CODE	DESCRIPTION	UOM
7666	Lastopress Cohesive Bandages 2.5cm x 2.75m Flesh	Each
7665	Lastopress Cohesive Bandages 5cm x 2.75m Flesh	Each
7664	Lastopress Cohesive Bandages 7.5cm x 2.75m Flesh	Each

### PORTABLE ICE MASSAGER

The Ice Up Portable Massager offers deep tissue relief for ligament, tendon and muscular injuries. Take it anywhere. Experience the benefits of direct, active ice massage:

- Effective treatments in 5-10 min. vs. 20 min. for passive ice packs
- Quick deep tissue penetration focused on target areas
- Effective therapy for ligament, tendon and muscular injuries
- Immediate ice massage increases treatment effectiveness, speeds recovery



CODE	DESCRIPTION	UOM
11080	Ice Up Portable Massager	Each

### WELCH ALLYN EXAM LIGHTS

#### GS300 General Exam Light

- A single LED that provides whiter, brighter light than halogen
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- A perfect combination of intense light output and cool operation all delivered in a compact design
- An affordable price and rugged performance
- Minimal heat dissipation
- Touchless on/off
- 5-year limited warranty



### USL SPORT WRAP 'N' ICE & REXIO ICE/HOT BAG

An easy and effective reusable way of applying ice and compression therapy to an injured area. The kit consists of a compression strap which fits over the top of an insulated rubber ice bag. Use with ice cubes, crushed ice or cold water.



CODE	DESCRIPTION	UOM
7458	USL Sport Wrap 'N' Ice	Each
7398	Rexio Care Ice/Hot Bag 11" 3000ml - Blue	Each

### EURO FACE, BATH & HAND TOWELS



CODE	DESCRIPTION	UOM
141232	EURO Face Cloth - White	Each
141231	EURO Hand Towel - White	Each
13985	EURO Bath Towel - White	Each

#### GS600 Minor Procedure Light

- Three LEDs that provide a whiter, brighter light than halogen
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- Wider head, designed specifically for demanding viewing conditions
- Uniform spot and superior brightness for minor procedures and exams in any healthcare environment
- Minimal heat dissipation
- Touchless controls with dimming features
- 5-year limited warranty

CODE	DESCRIPTION	UOM
4048960	Welch Allyn Light Mobile Stand Taller - GS300/600	Each
4044416	Welch Allyn Light GS300 Wall Mounted	Each
4044616	Welch Allyn Light GS600 Wall Mounted	Each

## ORDER BY EMAIL

### CUSTOMER SERVICE

customerservices@uslmedical.co.nz

## ORDER BY FREEPHONE

### AUCKLAND

0800 658 814

### CHRISTCHURCH

0800 804 546

## USL BIO-MEDICAL SERVICING

### AUCKLAND

0800 265 218

### EMAIL

uslservice@uslmedical.co.nz

## AUCKLAND OFFICE

494 Rosebank Road, Avondale

### AUCKLAND NUMBERS

Telephone: 09 829 0960

Facsimile: 09 829 0962

### CHRISTCHURCH NUMBERS

Telephone: 03 348 5260

Facsimile: 03 348 9830

[www.uslmedical.co.nz](http://www.uslmedical.co.nz)

**usl**  
universal  
specialist  
limited

Reception  
Showroom  
Customer  
Orders  
Distribution

**usl** aesthetics

SIGVARIS



# SIGVARIS COMPRESSION

Medical science meets high fashion and innovation



Compression therapy is a cornerstone in the treatment of venous disorders and diseases of the lymphatic system. Compression therapy is also used to prevent venous thrombosis.

SIGVARIS Medical Products, prescribed by your doctor and correctly fitted, may positively improve the health of your legs by reducing venous symptoms and improving the venous blood flow.

Wearing medical compression stockings regularly may improve your health and quality of life significantly.

