

EAR, EYES, NOSE & THROAT











WALL SETS & INTEGRATED WALL SETS

INTEGRATED WALL SYSTEMS

- Help reduce facility cost with an energy efficient product
- Standardise your facility with a sturdy, predrilled mounting board that mounts easily to any wall
- Choose either PanOptic™ Ophthalmoscope or the Coaxial Ophthalmoscope
- Choose either MacroView™ Otoscope or Diagnostic Otoscope
- Help improve workflow with contoured handles that automatically turn on/off
- Configure to meet your needs with the choice of the Wall Aneroid
- Optional KleenSpec® Disposable Dispenser
- Optional SureTemp® Plus 690 Electronic Thermometer



Coax	Ophthalmoscope		iExaminer	Wall Board		Wall Transformer	Aneroid Sphyg	Therm.		
	Basic	Plus		IWS	CSM/SPOT 44					
11720-L										
	118-2	118-3	238-2	238-3	11840-EX	411471	412260	77716	7670-01	01690-300
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1
1	1	1	1	1	1	1	1	1	1	1

777-PM3XXX-6
777-PM2XXX-6
777-SM2XXX-6
777-PM3WXX-6
777-PM2WXX-6
777-SM2WXX-6
777-PM3WCX-6
777-PM2WCX-6
777-PM3WAS-6
777-PM2WAS-6
777-SM2WAS-6
777-PM3WAX-6
777-PM2WAX-6
777-SM2WAX-6
777-PM3WXS-6
777-PM2WXS-6

GREEN SERIES™ 777 WALL TRANSFORMER

CONFIGURED TO EACH FACILITY'S NEEDS

- Help reduce facility cost with a more energy efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- Quickly and easily upgrade from the Welch Allyn® 767 Wall Transformer to the Green Series™ 777 Wall Transformer
- No on/off switch reduces contact, to help reduce the risk of cross-contamination



CODE	DESCRIPTION
4077716	W/A GS WALL TRANSFORMER WITH TWO HANDLES

* ADD THE OPHTHALMOSCOPE OR OTOSCOPE HEAD YOU WANT

OPHTHALMOSCOPE HEADS



CODE	DESCRIPTION
4011720-L	3.5V COAXIAL OPHTHALMOSCOPE LED
40118-2	PANOPTIC™ OPHTHALMOSCOPE - BASIC (20x viewing area)
40118-3	PANOPTIC™ PLUS (Sight align, iExaminer and connectivity options)

OTOSCOPE HEADS



CODE	DESCRIPTION
40250-2	LED OTOSCOPE
40238-2	MACROVIEW™ OTOSCOPE - BASIC
40238-3	MACROVIEW™ OTOSCOPE - PLUS - FOR IEXAMINER (iExaminer Smart Bracket purchased separately)

Optional extras

EAR SPECULUM DISPENSER



CODE	DESCRIPTION
4052400	EAR SPECULUM HOLDER/ DISPENSER

MOBILE STAND



CODE	DESCRIPTION
4400101	MOBILE STAND TO MOUNT THE WALL CHARGER ON

DESK CHARGER SET

UNIVERSAL DESK CHARGER

Pair your exams with our convenient, compact desk charging unit and exam tool sets. Drop-in desk/well charging Welch Allyn rechargeable power handles to help keep your otoscope, ophthalmoscope and other exam tools ready for use.



CODE	DESCRIPTION
40719DSK-6	W/A DESK CHARGER ONLY



CODE	DESCRIPTION
4971XXLDE-6	W/A DESK CHARGER WITH 2+ HANDLES

* ADD THE OPHTHALMOSCOPE OR OTOSCOPE HEAD YOU WANT

OPHTHALMOSCOPE HEADS



CODE	DESCRIPTION
4011720-L	3.5V COAXIAL OPHTHALMOSCOPE LED
40118-2	PANOPTIC™ OPHTHALMOSCOPE - BASIC (20x viewing area)
40118-3	PANOPTIC™ PLUS (Sight align, iExaminer and connectivity options)

OTOSCOPE HEADS



CODE	DESCRIPTION
40250-2	LED OTOSCOPE
40238-2	MACROVIEW™ OTOSCOPE - BASIC
40238-3	MACROVIEW™ OTOSCOPE - PLUS - FOR IEXAMINER (iExaminer Smart Bracket purchased separately)

Keep your diagnostic instruments together and easily portable, with our most popular handle and instruments in one easy-to-order set.

- Choice of the PanOptic™ Plus, Panoptic Basic or 3.5v Coaxial Ophthalmoscope Heads
- Choice of the Macroview™ Plus, Macroview Basic or LED Diagnostic Otoscope Heads
- Choice of the new Li-Ion Plus USB-C Handle, Standard Lithium Ion Handle or Convertible Silver handle uses C-cell batteries



CODE	DESCRIPTION
4071PM3LXE	PanOptic™ Plus, MacroView™ Plus, Li-Ion Plus USB-C, iExaminer™ - Hard Case- can be used to take images of ears and eyes.



CODE	DESCRIPTION
4071SM2LXU	LED Coaxial, MacroView™ Basic, Li-Ion Basic USB - Hard Case



CODE	DESCRIPTION
4071PM2LXU	PanOptic™ Basic, MacroView™ Basic, Li-Ion Basic USB - Hard Case



CODE	DESCRIPTION
4071SS2CX	LED Coaxial, LED Coaxial Ophthalmoscope, C-Cell Convertible Handle - Hard Case. (no batteries included uses C Cell)

PANOPTIC™ & PANOPTIC™ PLUS OPHTHALMOSCOPES offer a brighter, wider view.

Enhance eye exams with the PanOptic™ Plus, featuring QuickEye™ alignment technology and the option to save and share eye images. Both devices work on existing Welch Alllyn 3.5 V power sources.

USL Code	Description	Images	Examiner
40118-3	 <p>Welch Alllyn PanOptic™ Ophthalmoscope - Plus - For iExaminer (iExaminer Smart Bracket purchased separately) Quick Eye Alignment Technology</p>	 <p>Quick Eye™ Alignment - Two indicator lights help you direct the patient's gaze for better optic disc alignment</p>	
40118-2	<p>Welch Alllyn PanOptic™ Ophthalmoscope - Basic, (20x viewing area)</p>		 <p>Can't be used for taking pictures</p>

Specifications	PanOptic Plus Ophthalmoscope for iExaminer	PanOptic Basic Ophthalmoscope
Model Number	118-3	118-2
iExaminer ready	Yes	No
Quick Eye™ Alignment Technology	Yes	No
Warranty	5 Years	1 Year



3.5V STANDARD COAXIAL OPHTHALMOSCOPE HEAD

- Halogen HPX™ lamp provides 30% whiter, brighter light for true tissue colour and long-lasting performance
- Detect corneal abrasions and foreign bodies with cobalt blue filter
- LED bulb for brighter light



CODE	DESCRIPTION
4011720L	OPHTHALMOSCOPE LED COAXIAL-HEAD

3.5 V DIAGNOSTIC OTOSCOPES

Long-lasting LED illumination with fiber optics projects a cool light to the distal tip with no reflections or obstructions for superior views of the tympanic membrane.

Welch Allyn LED Otoscope




- LED technology for a brighter, whiter view
- Wide-angle viewing lens allows for instrumentation under magnification
- Sealed system for pneumatic otoscopy and accurate diagnosis of otitis media
- Works with all Welch Allyn 3.5V power sources

Welch Allyn MACROVIEW™ and MACROVIEW™ PLUS OTOSCOPES.


- Focus-free design
- Long-lasting LED illumination
- For brighter views, pair with our LumiView Clear Single-Use Ear Specula
- MacroView Plus otoscope models are digitally-capable, for capturing and sharing images of the ear

EXPERIENCE BIGGER, BRIGHTER VIEWS WITH OUR MACROVIEW™ AND MACROVIEW™ PLUS OTOSCOPES.

Enhance exams with the MacroView™ Plus Otoscope, and the option to capture and share ear images for consultation. Both devices work on existing Welch Allyn 3.5 V power sources.

USL Code	Description	Images	Examiner
40238-3	 <p>Welch Allyn MacroView™ Otoscope - Plus - For iExaminer (iExaminer Smart Bracket purchased separately)</p>	 <p>Pair the MacroView™ otoscope with Welch Allyn LumiView™ Clear Single-Use Ear Specula for brighter views.</p>	

Specifications	MacroView Plus Otoscope for iExaminer	MacroView Basic Otoscope
Model Number	238-3	238-2
iExaminer digital-ready	Yes	No
Warranty	5 Years	1 Year



40238-2	Welch Allyn MacroView™ Otoscope - Basic	  <p>Can't be used for taking pictures</p>
40250-2	Welch Allyn LED Otoscope	 

NEW LUMIVIEW™ CLEAR EAR SPECS

Features

- Disposable
- 4x brightness when compared to standard black specs
- Available in two sizes: Pediatric and Adult
- Compatible with most popular Welch Allyn® otoscopes
- Fit in 524 series spec dispensers that are currently in use
- Can be compatible with many competitive otoscopes

Benefits

- Disposable ear specs offer reduced risk of cross contamination
- Enhanced illumination and view of the ear canal and tympanic membrane (ear drum)
- Excellent user experience and feedback based on market research
- Upsell opportunity from our current disposable ear specs



CODE	DESCRIPTION
4052432CLR	W/A DISP CLEAR EAR SPEC 2.75mm (TUBE 34)
4052434CLR	W/A DISP CLEAR EAR SPEC 4.25mm (TUBE 34)



POCKET LED

SureColor LED Illumination

- **DURABILITY** - The most durable pocket set on the market based on an independent lab drop test study
- **LIGHTWEIGHT** - The lightest weight Welch Allyn® otoscope and ophthalmoscope ever
- **COLOURS** - The only set on the market available in a choice of four housing colours

CODE	DESCRIPTION
SETS IN SOFT CASE	
4092871BLU	POCKET SCOPE LED SET - BLUEBERRY
4092871BLK	POCKET SCOPE LED SET - ONYX
4092871PUR	POCKET SCOPE LED SET - PLUM
4092871WHT	POCKET SCOPE LED SET - VANILLA
OTOSCOPE ONLY	
4022880BLK	POCKET OTOSCOPE PLUS LED - ONYX IN CASE
4022880PUR	POCKET OTOSCOPE PLUS LED - PURPLE IN CASE
OPHTHALMOSCOPE ONLY	
4012880PUR	POCKET OPHTHALMOSCOPE + CASE - PURPLE
4012880BLK	POCKET OPHTHALMOSCOPE + CASE - BLACK
4012880WHT	POCKET OPHTHALMOSCOPE + CASE - WHITE
4012880BLU	POCKET OPHTHALMOSCOPE + CASE - BLUE

INDIVIDUAL HANDLES

INDIVIDUAL HANDLES

NB: If you want to purchase an individual Otoscope or Ophthalmoscope then add one of the below handles to one of the Otoscope or Ophthalmoscope heads on previous pages.



LITHIUM-ION HANDLE, CHARGING POD AND CHARGER

CODE	DESCRIPTION
4071906	BLACK LITHIUM-ION HANDLE + LITHIUM-ION CHARGING POD + CHARGER



NICKEL CADMIUM HANDLE

CODE	DESCRIPTION
4071020C + 4071036	SILVER NICKEL CADMIUM HANDLE PLUS CHARGER



CHARGER

LI-ION PLUS USB-C HANDLE

- Housing Design change
- Auto-Off feature (3mins)
- Low Battery Power Indicator
- New Design Rheostat - No Green Button
- USB-C Charging port
- Accepts existing 3.5 V heads
- Precise light control 5-stop setting
- Visually linear dimming
- 5.5 hours on-time when fully charged
- Can be charged using Desk Charger, Charging pod and USB-C charging cable. Fast Charge to 75% capacity in 20 minutes

CODE	DESCRIPTION
40719-3	LI-ION PLUS USB-C HANDLE



FIBRE OPTIC LARYNGOSCOPE BLADES

- Single-piece stainless steel construction ensures blade integrity
- Damaged light pathways can be repaired at a fraction of the cost of replacement blades
- Removable light pipes make these an affordable choice in fiber optic laryngoscope blades
- Upgrade your lightweight Fiber Optic Handles and Sets to SureColor LED technology. Code 406000-LED
 - Designed to enhance your procedures while offering superlative color rendition
 - Guaranteed for seven years

E MAC BLADE



CODE	DESCRIPTION
4069211	F/O - E MAC BLADE SIZE 1
4069212	F/O - E MAC BLADE SIZE 2
4069213	F/O - E MAC BLADE SIZE 3
4069214	F/O - E MAC BLADE SIZE 4

MACINTOSH BLADE



CODE	DESCRIPTION
4069061	F/O - MAC BLADE SIZE 1
4069062	F/O - MAC BLADE SIZE 2
4069063	F/O - MAC BLADE SIZE 3
4069064	F/O - MAC BLADE SIZE 4

MILLER STRAIGHT BLADE



CODE	DESCRIPTION
4068065	F/O - MILLER BLADE SIZE 00
4068060	F/O - MILLER BLADE SIZE 0
4068061	F/O - MILLER BLADE SIZE 1
4068062	F/O - MILLER BLADE SIZE 2
4068063	F/O - MILLER BLADE SIZE 3

FIBRE OPTIC HANDLES



CODE	DESCRIPTION
4060813	W/A FIBRE OPTIC HANDLE - C CELL
406014	W/A FIBRE OPTIC HANDLE - AA
4060815	W/A FIBRE OPTIC HANDLE - STUBBY



SPOT VISION SCREENER

The Spot Vision Screener is a handheld, portable device designed to help users quickly and easily detect vision issues on patients from 6 months of age through adult. Spot screens both eyes at once from a nonthreatening 3-foot distance. The touch-screen display allows for one-touch activation, simple management of patient data entry, and easy configuration for both vision screening and autorefractometry applications.

CODE	DESCRIPTION
40VS100S66	W/A SPOT VISION SCREENER

Available
to GP
customers
only

SCREENING AUDIOMETER

The Model 116 is a portable, manual screening audiometer designed for fast and efficient screening. It is used by primary care professionals around the world.

With the optional battery function and Amplivox Audiocups, the Model 116 can be used in a variety of normally challenging environments. The simple operation including the Threshold Retention function, allow the user to quickly store results during the test.

Results can then be stored on the internal memory, transferred to the PC or printer it includes:

- Audiometric headset
- Patient response switch
- Audiogram cards (50)
- Amplisuite Software - (allows you to enter recordings on the computer and print as a PDF)
- Operating manual
- Mains power adaptor
- Carrying case

INCLUDES
A PAIR OF
NOISE
REDUCING
HEADPHONES



Including



CODE 8507802

Optional items



CODE	DESCRIPTION
OMPT0	ECS Audiometer Thermal Printer
A108	Printer cable for audiometer
OM8029305	ECS Printer Paper - Sanibel Printer

OTOWAVE – TYMPANOMETER TWO OPTIONS

Available to GP customers only



OTOWAVE 102

The Otowave 102 is very intuitive Tympanometer to use with on-screen information displayed before, during and after the test. The compliance peak, compliance peak pressure, gradient, ear canal volume, The programmable reflex test can be configured for one frequency (1KHz) for the Otowave 102-1 - Tests can be sent to the computer using an Infrared reader and turned into a PDF

The Otowave 102 is also supplied with its own Data Management Software (ampliSuite) as well

COMES WITH

- Eartip selection box
- 4x test cavities
- Spare probe tips and gaskets
- 4 x AA batteries
- Carrying case
- AmpliSuite PC software (stand-alone)
- Noah impedance module

OTOWAVE 102-C

Lightweight and portable, the Otowave 102-C has the battery life to keep up with your whole working day.

With a charging and data transfer cradle, plus a user intuitive menu setup, the Otowave 102-C is perfect for any testing environment. Start the device, insert the probe and perform a tympanogram in only 3 seconds.

The Otowave 102-C is supplied with its own Data Management Software (ampliSuite).

COMES WITH

- Otowave 102-C
- Charging and Data Cradle
- Eartip selection box
- 4x test cavities
- Spare probe tips and gaskets
- Carry case
- AmpliSuite PC software (stand-alone)
- USB cable

Accessory for Otowave 102 to get tests into the computer




CODE	DESCRIPTION
OMT91	ECS Otowave IrDA PC Port

CODE	DESCRIPTION
OM001	Ecs Otowave 102 Tympanometer


Printer accessory for both units

CODE	DESCRIPTION
OMPT01	ECS Printer for Otowave
8029305	Paper for Sanibel Printer





OTOWAVE EARTIPS
(single patient use)



CODE	DESCRIPTION
OMIA011	OTOWAVE EARTIP 11mm RED
OMIA012	OTOWAVE EARTIP 12mm BLUE
OMIA013	OTOWAVE EARTIP 13mm GREEN
OMIA014	OTOWAVE EARTIP 14mm YELLOW
OMIA015	OTOWAVE EARTIP 15mm BLUE
OMIA019	OTOWAVE EARTIP 19mm YELLOW

Example of PDF

COMPARE USER FEATURES

	Otowave 102	Otowave 102-C
Totally portable	✓	✓
Cradle option Nil	Nil	✓
PC connectivity (USB)	Nil	✓
PC connectivity (IrDA)	✓	Nil
Thermal printer output	✓	✓
Internal memory	✓	✓
Battery power function	✓ AA Batteries	✓ Built-in rechargeable battery
Mains power function (adaptor supplied)	Nil	✓ via cradle
Carry case supplied	✓	✓

BECAUSE WE CARE.



USL CONTACT

DEBBIE HOLLIS
PRODUCT MANAGER
- WELCH ALLYN
debbie@uslmedical.co.nz
027 486 9589

MOZINA NISHA
EQUIPMENT & MARKETING
CO-ORDINATOR
mozina.nisha@uslmedical.co.nz

0800 658 814 | www.uslequipment.co.nz

 **usl** equipment