

# THERMOMETRY



## PRO 6000 EAR THERMOMETER

The Welch Allyn® Braun ThermoScan® PRO 6000 ear thermometer featuring advanced measurement technology, has the speed and simplicity you need to capture accurate and reliable temperature readings in a comfortable way for patients.

The Braun PRO 6000 advantage lies in our innovative PerfectTemp™ technology, which overcomes the potential for low readings compared to core, by adjusting for factors that impact measurement accuracy.

### FEATURES

- Innovative PerfectTemp™ technology adjusts for variability in probe placement
- ExacTemp™ technology detects stability of the probe during measurement
- Designed to take quick and accurate readings, to improve patient satisfaction
- 60 second pulse timer assists with manual measurement of pulse rate, respiration
- Electronic security features help prevent theft, loss
- Small or large cradles for probe cover storage that pop-up for easy attachment
- 6 foot or 9 foot security tether options keep device attached to cradle



PRO6000 WITH CRADLE



2.7M TETHER



1.8M TETHER



PR6000 WITH LARGE CRADLE



PROBE COVERS BX200

CODE	DESCRIPTION
406000200	THERMOSCAN BRAUN PRO6000 WITH CRADLE
406000300	THERMOSCAN BRAUN PRO6000 WITH LARGE CRADLE

CODE	DESCRIPTION
40106204	2.7M TETHER
40106201	1.8M TETHER
4060005	THERMOSCAN PROBE COVERS BX200



## SURETEMP® PLUS 690

Interchangeable, removable probe wells reduce the risk of cross-contamination

- Large LCD displays temperatures in Fahrenheit or Celsius
- Last temperature recall
- Convenient storage for 25 probe covers
- Takes approximately 6,000 readings on 3 AA batteries



## SURETEMP® PLUS 692

- Ergonomic design - comfortable to hold, easy to use, and very durable
- Interchangeable, removable probe wells reduce the risk of cross-contamination
- Large LCD display with location identification capability
- Includes a 60-second pulse timer for counting pulses and respirations
- Mechanical or electronic security options to help prevent theft
- Takes approximately 5,000 readings on three AA batteries
- Wall mount standard with the SureTemp® Plus 692 model