




**UVA** equipment

# EQUIPMENT PRODUCT GUIDE





## OSCIMED CAST SAWS AND EXTRACTION UNITS

The Oscimed 2000 is a very low noise level cast saw with patented anti-noise technology. The Oscimed 2000 consists of a brushless, fanless 100W-24V motor.

It has been ergonomically designed with a modern shape for single handed use. The saw has an instant blade change ability with a patented non-slip lock nut.

Oscimed technology includes saw blades with a noise-suppressing central core.

This vacuum cleaner possesses four levels of filtration (including two filters to HEPA standards) and filters dust down to 3 microns. It has a tough, powerful construction with a support for the saw and vacuum hose. The unit comes complete with an electricity supply socket for the saw. Available in 1-speed version (1000W).

- Sound level: 62dB in speed 1
- Suction: 2200mm water gauge, 50 litres/sec
- Length of vacuum hose: 3 metres
- Capacity of stainless steel tank: 20 litres
- Weight: 12 kg



### OSCIMED CAST SAW AND EXTRACTION VACUUM UNIT

CODE	DESCRIPTION
20101	Oscimed Cast Saw
20102	Oscimed Extraction Vacuum Unit
20102B	Oscimed Blade For Synthetic Blade 65mm
20102C	Oscimed Blade for POP 65mm
20102D	Oscimed Blade - Titanium 65mm
20102A	Oscimed Dust Bags 10/box
13413	Oscimed Cast Spreader



## PROUDLY NZ OWNED

USL has been built by creating seamless relationships. Since 1984 we've built enduring partnerships with customers, suppliers and our community. We set out to create the very best environment for our team so that they might deliver their unique skills to our markets.

### WHAT YOU CAN EXPECT FROM US:

- We offer our customers dependability and experience they can rely on and trust.
- We deliver under IFOTIS: In-Full, On-Time, In-Specification.
- We operate a management reporting system that provides information on the frequency, recency, amount and type of purchases for all customers, and thereby enhancing their business purchasing efficiency.
- Normal guarantees apply and are rigorously honoured, and we guarantee to supply our key clients with critical stock items under IFOTIS.
- Our product range is comprehensive. It meets our Healthcare customers' needs, so they can confidently have a one supplier relationship for all medical products.
- Private Market
- General Practitioners
- Nurse Practitioners
- Private Hospitals (Medical, Surgical)
- Rest Homes
- Accident and Emergency Clinics
- Emergency Services
- Occupational Health Services
- Community Trusts
- Physiotherapists and Sports Medicine Clinics
- Sports Clubs and Organising Bodies
- Veterinary Services
- First Aid Suppliers
- Medical Laboratories
- Public Hospitals
- Specialist Clinics
- Community and Public Health Services
- Pharmacy Retail
- Industrial and Safety
- Dermatologists and Skin Clinics

### OUR COMPANY AT A GLANCE:

- USL is a limited liability company, proudly New Zealand owned
- Our annual sales exceed \$NZ 60 Million
- Our head office and main distribution centre are based in Auckland
- The South Island is serviced from our Christchurch distribution centre
- We service 3,750 customers at all locations throughout New Zealand including all 20 of the National District Health Board



Order By Freephone:  
Auckland: 0800 658 814  
Christchurch: 0800 804 546

Order By Freefax:  
Auckland: 0800 830 660  
Christchurch: 0800 804 547

Auckland Head Office  
494 Rosebank Road, Avondale, Auckland 1026  
Email: [customerservices@uslmedical.co.nz](mailto:customerservices@uslmedical.co.nz)